

HSW ECO-MATIC® TWIN

0,5 ml / 1 ml / 2 ml

Automatic self-filling syringe with two separate barrels
Jeringa automática con dos cilindros separados

- Premium-quality plastic syringe with durable plastic barrels
 - Ergonomic, light weight design for ease of use
 - Simultaneous vaccination of two products independently
 - Bottle mount and/or tubing connection
 - Easy adjustable needle distance
 - Easy adjustable barrel position
 - Various fixed dose barrel combinations available
-
- Jeringa de plástico de primera calidad con cilindros de plástico de gran rigidez
 - Diseño ergonómico, manejo fácil a pesar de un peso ligero
 - Administración de dos medicamentos / vacunas
 - Adaptador de botella y/o adaptador de tubo son posibles
 - Distancia de la aguja es ajustable
 - Posición del cilindro es ajustable
 - Dosificación fija, disponible en diversas combinaciones



Included in delivery

HSW ECO-MATIC® TWIN complete:

1 x Syringe
2 x Universal bottle collars Ø 20 mm
2 x Universal bottle collars Ø 30 mm
2 x Tubing collars (Ø 4-6 mm)
2 x PVC tubing (Ø 4 mm)

Alcance del suministro

HSW ECO-MATIC® TWIN completo:

1 x Jeringa
2 x Adaptador universal Ø 20 mm
2 x Adaptador universal Ø 30 mm
2 x Adaptador de tubo (Ø 4-6 mm)
2 x Dos tubos PVC (Ø 4 mm)



Version
Modelo

Vol/ml

Item no.
No. de artículo

Luer Lock	0,5 ml + 0,5 ml 0,5 ml + 2 ml 1 ml + 1 ml 2 ml + 2 ml 1 ml + 2 ml	8300037936 8300031731 8300037673 8300037669 8300037671
-----------	---	--



The needle distance can be adapted to the age, size and species of the animal.
There are seven different distances:
20 mm - 21,1 mm - 24 mm - 28 mm - 32 mm - 34,9 mm - 36 mm

La distancia entre las agujas puede adaptarse a la edad, la tamañ a y a la especie del animal. Hay 7 distancias diferentes:
20 mm - 21,1 mm - 24 mm - 28 mm - 32 mm - 34,9 mm - 36 mm

